LS3 3mm Uni-Stable Shutter Specifications


Features

- The UNIBLITZ® LS3 series is a 3mm clear aperture especially well suited for laser use.
- Designed for low level chopping, pulse gating, selection, and modulation to 200Hz.
- Precision exposure control in the areas of holography and laser switching are additional applications that benefit from the precise, repeatable characteristics of the LS shutter series.
- Three aperture sizes available: 2mm, 3mm, and 6mm.
- Laser energy ratings up to 5W/mm2 with "Z" or "ZM" shutter blade coating options.
- Exposure repetition rates from DC 200Hz.
- Electronic Synchronization System available.
- Available in a normally-open configuration.
- Available in an un-housed configuration through special order.



Fig. 1 LS3 3mm Uni-Stable Shutter

Timing



¹Typical timing values (msec.) using UNIBLITZ® drive equipment and measured with UNIBLITZ® shutters equipped with standard black Teflon® coated shutter blades.

	LS3	Time (msec.) ¹
O-A:	Delay time on opening after current is applied	1.0
A-C:	Transfer time on opening	0.5
O-C:	Total opening time	1.5
C-E:	Min. dwell time with min. input pulse	0.8
B-F:	Min. equivalent exp. time	1.5
E-G:	Transfer time on closing	0.6
A-G:	Total window time	1.9
MET:	Min. exposure time	2.0
TEP:	Typical exposure pulse	>2.3

⁸⁰³ Linden Avenue • Rochester, NY 14625 • 585-385-5930 tel • 800-828-6972 toll-free • 585-385-6004 fax • www.uniblitz.com



Electrical

COIL RESISTANCE	PULSE VOLTAGE TO OPEN	HOLD VOLTAGE ¹ (NOMINAL)
48 OHMS	+65 VDC	+10 VDC

¹ Voltage level required across actuator coil when being held in the open position.

The Electronic Synchronization System provides a feed-back signal (through the driver utilized) after the shutter transfers to the open state. The system incorporates an infrared emitting diode, an infrared sensitive detecting transistor, and an interrupting vane. The vane is attached to the shutter so as to block the light path between the emitter and detector in the closed position. When the shutter transfers to the 80% open position, the vane is removed from the infrared light path, allowing the emitter to switch the detector to the active state. No connection to the designated synchronization pins when no electronic sync. is selected.




Fig. 3

Mechanical


SERIES	WEIGHT UNCASED	WEIGHT CASED	OPERATING TEMP. (DEGREES)	MAX. OPENING BOUNCE	MAX. CLOSING BOUNCE	MAX. FREQUENCY OF OPERATION ²	NUMBER OF SHUTTER BLADES
LS3	N/A	7.41 oz. (.21 Kg)	0-80 ° C	15%	5%	50 Hz / 200 Hz	2

² (CONT/BURST) CONTinuous frequency rating specified at shutter's minimum exposure pulse. BURST frequency rating specified for (4) four seconds maximum with (1) one minute minimum between bursts. Frequency measurements are taken in free air, 25°C ambient, actuator coil equipped with heat sink. For additional information on maximum sustained frequencies obtainable, please contact one of our technical representatives.

803 Linden Avenue • Rochester, NY 14625 • 585-385-5930 tel • 800-828-6972 toll-free • 585-385-6004 fax • www.uniblitz.com



Housing



- The LS #2 case style allows a number of mounting configurations.
- A ½-20threaded hole is provided for post mounting.
- The 1.625 in x 40TPlexternal thread, rear side, and the six 2-56 threaded holes, front side, can be interfaced directly into your application or fitted with a variety of specific mounting options
- See "MICROSCOPE & VIDEO MOUNTING SYSTEMS" data sheet for additional information.
- The unit terminates with a 7-pin male connector.
- Drawings of the device in its normally open and un-housed configuration are available under the 'Downloads' tab on our website.

Fig. 4

Housing/Connector Layout




Fig. 5 illustrates the 7-pin connector and 1/4-20 threaded hole layout for the LS2/3/6 series #2 case style.

803 Linden Avenue • Rochester, NY 14625 • 585-385-5930 tel • 800-828-6972 toll-free • 585-385-6004 fax • www.uniblitz.com



Optical Blade Finish

SHUTTER	ULTRAVIOLET .34µm (microns)		VISIBLE .475µm (microns)		INFRARED .75- 10.6µm (microns)	
SERIES	(Z) AlSiO	(ZM) AIMgF ₂	(Z) AlSiO	(ZM) AIMgF ₂	(Z) AlSiO	(ZM) AIMgF ₂
LS3	N/A	2.5 W/mm ²	5 W/mm ²	2.5 W/mm ²	2.5 W/mm ²	2.5 W/mm ²

For reflectance graph, please visit http://www.uniblitz.com/optical-shutters-comparison-chart.aspx

Blade Samples are available upon request.

Product Options

Part Number: LS3S @ @ @ @ @ — @ — @ — @

Example Part Number: LS3S2T0R3-EC-NL-100

Aperture Size\Type: LS3S: 3mm (Normally closed) LS3E: 3mm (For use with D880C Driver) LSR3S: 3mm (Normally open) LSR3E: 3mm (Normally open, for use with D880C Driver)	2 Housing: 5 2: #2 Housing 4 Electronic Sync: 0: Omit Sync. 1: Electronic Sync. Included	3 Blade Finish: 3 T: Teflon® Coated S.S. Blades ZM: AlMgF2 Coated BeCu Blades Z: AlSiO Coated BeCu Blades
 High Temp Mod: R3: High Temp. Modification² Leave blank if not required 	 Encapsulated Coil EC: Encapsulated Coil Included² (For use with vacuum) Leave blank if not required 	RoHS Compliant Version: NL: RoHS Compliant ² • Leave blank if not required

- 8 Mounting Options (by type, #2 Housing Required):
 - Leave blank if not required
 - 21: Zeiss Axiovert Type 22: Nikon SBX Type 23: Olympus BH/IMT Type
 - 24: Olympus BX/IX Type 26: Leica Type 27: Nikon Type 28: Olympus IX Transmitted Type
 - 29: Nikon TE Type 30: Leica DM/DMIR/DMIRB Type 31: Nikon Confocal Type
 - 32: Nikon 80i Type 100: Mounting Ring 105: C-Mount Adapter (Male)
 - 106: C-Mount Adapter (Female) 110: T-Mount Adapter

For information regarding applicable intellectual property, please visit www.uniblitz.com/company-info/patents.

803 Linden Avenue • Rochester, NY 14625 • 585-385-5930 tel • 800-828-6972 toll-free • 585-385-6004 fax • www.uniblitz.com

¹ Input side only, Teflon® coating is on opposite side. Intended to protect the shutter blade surface, light source must be input to the reflective side only.

² Please visit our website for more information regarding this option.

³ Other blade options are available through special order.

⁴ If #2 Housing is selected, modification to housing is required.

⁵ An un-housed configuration is available through special order.